

Instructions for Conducting NMR-Based Kinetics Measurements

Laura Steffens, 5/2010

Slight Editing by P.F. Flynn

- I. Setting Up the Experiment
 - A. Optimize lock and shim at desired temperature with representative sample in NMR (ie. proper concentrations for actual NMR run)
 - B. Load one-dimensional NMR experiment, typically ^1H
 - C. Set spectral parameters (nt, np, sw, ss, etc.), appropriately
 - D. Set **pad** equal to zero (pad=0)
 - E. Type **time** <return> to obtain time for one experiment (e.g. ~347 s for one pulse sequence for a 20 mM sample)
 - F. Decide the overall length of the kinetics run (e.g. 4 hrs (2400s)) and at what increments you would like to record spectra (e.g. every 10 mins (600s)) with 24 increments (2400s/600s)
 - G. Determine preacquisition delay (pad) according to the length of one pulse sequence and increment size (e.g. 600s-347s=253s=pad)

- II. To create array
 - A. Type **array** <return>
 - B. Instrument prompts *param to array*, type **pad** <return>
 - C. Instrument prompts *# steps in array*, type **24** <return>
 - D. Instrument prompts *starting value*, type **253** <return>
 - E. Instrument prompts *array increment*, type **0** <return>
 - F. Need to adjust first delay
 - i. Type **pad[1]=0** <return>
 - G. Type **time** <return> to get (approximate) total time for arrayed experiment
 - H. Type **ga** <return> to begin experiment

- III. Processing Spectra
 - A. Type **ai** to set display in absolute intensity mode. for coarse baseline correction
 - B. Type **cdc dc** for a coarse baseline correction (see application note on baseline correction for a better way to do this – if necessary)
 - C. Viewing Spectra
 - i. Type **ds(1)** to load first spectrum
 - ii. Type **vsadj** to scale first spectrum (note that the integral intensity does not depend on vs, but it is useful to review the data set)
 - iii. Type **dssh** to see comparison of spectra (inspect relative intensities)
 - D. Printing
 - i. To print one spectrum, type normal print commands
 - ii. To print all spectra in a horizontal fashion, type **pl('all') page**

E. Saving

- i. To save fids, use normal save commands
- ii. To save display array of spectra, type **pl('all') page('filename.ps')**
- iii. To move saved array to archive
 1. Type **cp *filename* /archive/*your archive name*** <return>